

## STAINING F-ACTIN WITH PHALLOIDIN

### Materials and Reagents:

#### *S mix*

	[stock]	Volume added (µL)	[final]
Na-Phosphate buffer pH 7.5	0.8M	250	0.2mM
MgCl <sub>2</sub>	1M	1	1mM
SDS	1%	4	0.004%
dH <sub>2</sub> O	-	745	-

Na-phosphate: 8.1 ml 1M Na<sub>2</sub>HPO<sub>4</sub>, 1.9 ml 1M NaH<sub>2</sub>PO<sub>4</sub>, 2.5 ml dH<sub>2</sub>O

PBBT: PBS, 0.5% BSA, 0.5% Tween-20

Mounting solution: 90% glycerol, 10% PBS, 1 mg/ml phenylenediamine  
Rhodamine-phalloidin (Molecular Probes)

### Procedure:

- 1) Wash worms off the plate in water, PBS, or M9.
- 2) Wash 2 times to remove excess bacteria (spin down at 3000rpm X 1min).
- 3) Remove as much supernatant as possible.
- 4) Add 10µL worms to an eppendorf tube.
- 5) Freeze tube in liquid nitrogen.
- 6) Lyophilize worms by spinning in speed vac (~5 minutes).
- 7) Add 3-4 drops ice cold acetone and put in freezer 3 min.
- 8) Remove acetone with pipetman (or aspirate carefully) and speedvac to remove residual acetone.
- 9) Add 2U rhodamine-phalloidin to an eppendorf tube and speedvac to evaporate methanol. Resuspend in 20µL S mix.
- 10) Resuspend dry worms in S mix/phalloidin.
- 11) Stain in the dark at RT for at least 30 min.
- 12) Wash 2X in 1 mL PBBT.
- 13) Resuspend in ~20 µL PBS.
- 14) For viewing, mix equal volumes of worms and mounting medium.
- 15) Mount on an agar pad and examine under the microscope.

**Tips/Troubleshooting:**

- 1) Worms often remain stuck together after lyophilizing. Mix by pipeting up and down with a pipetman, but be sure to cut the tip to give a wider opening so that worms will not be sheared
- 2) Mounting medium is not strictly necessary unless you plan to store the slides for some time (days). Worms can be mounted in PBS on an agar pad.
- 3) Molecular Probes sells phalloidin conjugated to different fluorophores depending on what you need to co-stain with

**References**

- 1) This protocol is adapted from Michael Koelle's (<http://cobweb.dartmouth.edu/~ambros/worms/index.html>)

**Submitted by:**

Jim Morley 5/03